

V-Chromer Mini

Manual stainer for microscopic samples



Manual stainer with 6 staining containers



Convenient digital timer
(count-up and -down) included

Main features

- Use reagents of various manufacturers
- Resistant to chemical substances and stains
- Only 80 ml of reagent is used per container
- Reagents are protected from evaporation due to caps on containers and sample holders

Staining in only 3 simple steps



1

Place slides
in the holder



2

Move the holder from one
container to the next
according to the staining
protocol



3

Get standardized
stained slides

Pre-classification has to be validated by the specialist reviewing the galleries of cells. We reserve the right to change specification without notice.

